



Information about data generated by the device (according to the Data Act)

Product category: Treadmill
Product name: Zipro Chill

Data generated by the device

Data is generated to enable the user to monitor training parameters.

Parameter	Data generation method	Data storage location	Data retention period	Data format	Data volume	User access to data	Data sharing	Data export capability	Possibility to delete saved data
Time	The device collects data continuously during training.	Local computer memory (EEPROM)*	Locally on the device – until a new training session begins or the device is switched off.	FIT (Flexible & Interoperable Data Transfer)	One packet: approx. 10–20 bytes. 1 hour of training: approx. 18–36 kB. A session lasting approx. 3 hours =< 100 kB.	Via the computer, only until a new training session begins or the device is switched off.	Integration with applications (list of supported applications below).	No	-
Distance									
Speed									
Calories									
Body Fat	The device calculates data based on parameters entered by the user: gender, age, weight.	Local computer memory (EEPROM)*	Locally on the device – until a new training session begins or the device is switched off.	Integer or floating-point number, stored in binary format.	2 bytes	Via the computer, only until a new training session begins or the device is switched off.	No	No	-

* Data is saved locally, anonymously, and is not assigned to a specific user.

Applications

The user has full control over the use of applications. Using the treadmill does not require application installation. Using applications requires acceptance of the privacy policy.

Application name	Data download method	Data storage location	Data retention period	Data format	Data volume	User access to data	Data sharing	Data export capability	Ability to delete saved data
FitShow	Training and telemetry data (e.g. speed, cadence, heart rate) collected from equipment via BLE or sensors.	In the cloud or on servers Fitshow (in some versions data is stored in the People's Republic of China).	As long as necessary for the purposes for which they were collected or according to legal regulations.	FIT (Flexible & Interoperable Data Transfer)	-	The user has access to training session history; they can view their data in the application.	Possible disclosure in accordance with legal requirements. External service providers (e.g. analytics, notifications), synchronization	No	Yes — the user can delete their account / activities; data is deleted or made anonymous.

							with external accounts if the user enables such an option.		
Kinomap	Training and telemetry data (e.g. speed, cadence, heart rate) collected from equipment via BLE or sensors.	Data stored on Kinomap servers.	Up to 10 years, for some training data / user accounts.	Second-by-second telemetry data	-	The user has access to their activity history, can set privacy settings (public/private), and export their files (e.g. GPX).	Data may be shared: - with external platforms upon request or by the user (synchronization) - with technical service providers - for marketing/analytical purposes - in case of legal requirements.	Yes — GPX track export, activity sharing capability.	Yes — the user can delete their account / activities; data is deleted or made anonymous.
Zwift	Data provided by the user, collected from sensors and activities, or calculated by the application.	Data is stored on Zwift servers locally on the device.	Stored while the account is active and in accordance with legal requirements; no single established period.	Mainly FIT files, logs, and training files.	Variable; depends on the number and duration of activities.	Users can view profile data and activity history in the application.	Shared with other users, service providers, and connected platforms, depending on settings.	Activities can be exported as FIT files; users can also request a copy of their data.	Users can delete individual activities or their entire account.